

Giardia lamblia
10E7

Certificate of Analysis

Batch number: 480 - 254

Production Data

Volume: 2.45ml
Viability: Viable
Catalogue Number: G10E7-V
BTF Production Date: 22 April 2014

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 10 April 2014
Date Excreted: 10 April 2014
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Manager:



.....
Sue Warner

Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

