

Cryptosporidium parvum
10E7

Certificate of Analysis

Batch number: **482 - 258**

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: C10E7
BTF Production Date: 19 May 2014

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 8 April 2014
Date Excreted: 8 April 2014
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Manager:



.....
Sue Warner

Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

