

**Giardia lamblia**  
10E7

## Certificate of Analysis

**Batch number: 485 - 256**

### Production Data

**Volume:** 2.5ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** G10E7  
**BTF Production Date:** 19 May 2014

### Stock specifics

**Strain:** H3  
**Source:** Gerbil  
**Date purified:** 7 May 2014  
**Date Excreted:** 7 May 2014  
**Purification Method:** Sucrose and Percoll density gradient centrifugation

### Quality Manager:



.....  
Sue Warner

### Water Industry Enquires:

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

