

**Cryptosporidium parvum**  
10E7

Certificate of Analysis

Batch number: **508 - 274**

**Production Data**

Volume: 1ml  
Viability: Viable  
Catalogue Number: C10E7-V  
BTF Production Date: 12 November 2014

**Stock specifics**

Strain: IOWA  
Source: Bovine  
Date purified: 18 October 2014  
Date Excreted: 18 October 2014  
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

**Quality Officer:**



.....  
Linh Le

**Water Industry Enquires:**

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

