

**Cryptosporidium parvum**  
10E7

## Certificate of Analysis

**Batch number: 509 - 275**

### Production Data

**Volume:** 1ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** C10E7  
**BTF Production Date:** 1 December 2014

### Stock specifics

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 14 November 2014  
**Date Excreted:** 14 November 2014  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

### Quality Officer:



.....  
Linh Le

### Water Industry Enquires:

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

