

**Cryptosporidium parvum**  
10E7

## Certificate of Analysis

**Batch number: 518 - 286**

### Production Data

**Volume:** 1ml  
**Viability:** Viable  
**Catalogue Number:** C10E7-V  
**BTF Production Date:** 2 March 2015

### Stock specifics

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 14 December 2014  
**Date Excreted:** 14 December 2014  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

### Quality Officer:



.....  
Linh Le

### Water Industry Enquires:

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

