

**Cryptosporidium parvum**  
10E7

## Certificate of Analysis

**Batch number: 548 - 310**


### Production Data

**Volume:** 1ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** C10E7  
**BTF Production Date:** 10 August 2015

### Stock specifics

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 3 July 2015  
**Date Excreted:** 3 July 2015  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

### Quality Assurance Officer:



.....  
Fiona Pollock

### Water Industry Enquires:

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

