

Cryptosporidium parvum
10E7

Certificate of Analysis

Batch number: 549 - 315

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: C10E7
BTF Production Date: 7 September 2015

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 21 July 2015
Date Excreted: 21 July 2015
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Officer:



.....
Fiona Pollock

Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

