

Giardia lamblia

10E6

Certificate of Analysis

Batch number: **554 - 319**

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: G10E6
BTF Production Date: 6 October 2015

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 23 September 2015
Date Excreted: 23 September 2015
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Officer:



.....
Linh Le

Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

