

Cryptosporidium parvum
10E7

Certificate of Analysis

Batch number: 561 - 324

Production Data

Volume: 1ml
Viability: Viable
Catalogue Number: C10E7-V
BTF Production Date: 30 November 2015

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 20 October 2015
Date Excreted: 20 October 2015
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Officer:



.....
Fiona Pollock

Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

