

Cryptosporidium parvum
10E7

Certificate of Analysis

Batch number: 600 - 354

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: C10E7
BTF Production Date: 17 August 2016

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 30 June 2016
Date Excreted: 30 June 2016
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Officer:



.....
Fiona Pollock

Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

