

**Giardia lamblia**

10E6

Certificate of Analysis

Batch number: **607 - 361**

**Production Data**

Volume: 1ml  
Viability: Inactivated - Gamma Irradiation  
Catalogue Number: G10E6  
BTF Production Date: 14 September 2016

**Stock specifics**

Strain: H3  
Source: Gerbil  
Date purified: 31 August 2016  
Date Excreted: 31 August 2016  
Purification Method: Sucrose and Percoll density gradient centrifugation

**Quality Assurance Officer:**



.....  
Fiona Pollock

**Water Industry Enquires:**

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

