

**Cryptosporidium parvum**  
10E7

## Certificate of Analysis

**Batch number: 624 - 374**

### Production Data

**Volume:** 1ml  
**Viability:** Viable  
**Catalogue Number:** C10E7-V  
**BTF Production Date:** 30 January 2017

### Stock specifics

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 6 January 2017  
**Date Excreted:** 6 January 2017  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

### Quality Assurance Officer:



.....  
Fiona Pollock

### Water Industry Enquires:

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

