

**Cryptosporidium parvum**  
10E9

## Certificate of Analysis

**Batch number: 690 - 412**

### Production Data

**Volume:** 2ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** C10E9  
**BTF Production Date:** 29 March 2018

### Stock specifics

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 21 March 2018  
**Date Excreted:** 21 March 2018  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

### Quality Officer:

  
.....

Linh Le

### Water Industry Enquires:

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

