

Giardia lamblia

10E6

Certificate of Analysis

Batch number: **697 - 418**

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: G10E6
BTF Production Date: 23 May 2018

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 11 May 2018
Date Excreted: 11 May 2018
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Officer



.....
Madhuli Patil

Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

