

**Giardia lamblia**

10E6

Certificate of Analysis

Batch number: **701 - 421**

**Production Data**

Volume: 1ml  
Viability: Inactivated - Gamma Irradiation  
Catalogue Number: G10E6  
BTF Production Date: 19 June 2018

**Stock specifics**

Strain: H3  
Source: Gerbil  
Date purified: 5 June 2018  
Date Excreted: 5 June 2018  
Purification Method: Sucrose and Percoll density gradient centrifugation

**Quality Assurance Officer**



.....  
Madhuli Patil

**Water Industry Enquires:**

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

