

***Cryptosporidium parvum***  
**10E7**

Product Information Sheet

**Batch number: 722 - 442**

**Production Data**

**Volume:** 1ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** C10E7  
**BTF Production Date:** 15 October 2018

**Stock specifics**

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 29 September 2018  
**Date Excreted:** 29 September 2018  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Quality Manager**



.....  
**Peter Verwoerd**  
19 October 2018



**Water Industry Enquires:**

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8221 9589  
[www.btfbio.com](http://www.btfbio.com)



Issue: 1 - August 2018