

Cryptosporidium parvum
10E7

Product Information Sheet

Batch number: 725 - 444

Production Data

Volume: 1ml
Viability: Viable
Catalogue Number: C10E7-V
BTF Production Date: 30 October 2018

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 27 September 2018
Date Excreted: 27 September 2018
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Manager



.....
Peter Verwoerd
2 November 2018



Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8221 9589
www.btfbio.com

