

Cryptosporidium parvum
10E8

Product Information Sheet



Batch number: 730 - 447

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: C10E8
BTF Production Date: 5 December 2018

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 24 October 2018
Date Excreted: 24 October 2018
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Manager

.....
Peter Verwoerd
13 December 2018



Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8221 9589
www.btfbio.com



Issue: 1 - August 2018