

# ***Giardia lamblia***

# Product Information Sheet

10E6

**Batch number: 732 - 450**

## **Production Data**

**Volume:** 1ml  
**Viability:** Inactivated - Gamma Irradiation  
**Catalogue Number:** G10E6  
**BTF Production Date:** 5 December 2018

## **Stock specifics**

**Strain:** H3  
**Source:** Gerbil  
**Date purified:** 23 November 2018  
**Date Excreted:** 23 November 2018  
**Purification Method:** Sucrose and Percoll density gradient centrifugation

## **Quality Manager**



.....  
**Peter Verwoerd**  
13 December 2018



## **Water Industry Enquires:**

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8221 9589  
[www.btfbio.com](http://www.btfbio.com)



Issue: 1 - August 2018