

**Cryptosporidium parvum**  
10E8

Product Information Sheet

Batch number: **733 - 454**

**Production Data**

**Volume:** 1ml  
**Viability:** Viable  
**Catalogue Number:** C10E8-V  
**BTF Production Date:** 3 January 2019

**Stock specifics**

**Strain:** Iowa  
**Source:** Bovine  
**Date purified:** 18 November 2018  
**Date Excreted:** 18 November 2018  
**Purification Method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Quality Assurance Specialist**



.....  
**Linh Le**  
4 January 2019



**Water Industry Enquiries:**

BTF Pty Ltd.  
Sydney Australia  
Tel. 61 (0) 2 8877 9150  
Fax. 61 (0) 2 8221 9589  
[www.btfbio.com](http://www.btfbio.com)



Issue: 1 - August 2018