

Giardia lamblia
10E6

Product Information Sheet

Batch number: 735 - 451

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: G10E6
BTF Production Date: 2 January 2019

Stock specifics

Strain: H3
Source: Gerbil
Date purified: 14 December 2018
Date Excreted: 14 December 2018
Purification Method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Specialist



.....
Linh Le
14 January 2019



Water Industry Enquiries:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8221 9589
www.btfbio.com



Issue: 1 - August 2018