

Cryptosporidium parvum
10E7

Product Information Sheet

Batch number: 740 - 455

Production Data

Volume: 1ml
Viability: Inactivated - Gamma Irradiation
Catalogue Number: C10E7
BTF Production Date: 29 January 2019

Stock specifics

Strain: Iowa
Source: Bovine
Date purified: 12 January 2019
Date Excreted: 12 January 2019
Purification Method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Specialist



.....
Linh Le
6 February 2019



Water Industry Enquiries:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8221 9589
www.btfbio.com



Issue: 1 - August 2018