

**Batch Number: B 196**

**General**

**Safety:** This product is non-hazardous.  
**Storage:** 2-8°C. Do not freeze.  
**Catalogue Number:** CSCG100  
**Preparation Date:** 13 November 2014  
**Expiration Date:** 13 March 2015  
**Volume:** 1.098 ml +/- 3 ul  
**Suspension media:** Buffered saline solution.  
**Sterilisation method:** Gamma irradiation.

<b>Counts</b> <small>(Method Ref: CG-011)</small>		
	Mean	SD
<b>Cryptosporidium count:</b>	<b>100</b>	<b>1.5</b>
	Mean	SD
<b>Giardia count:</b>	<b>99</b>	<b>1.4</b>
<b>DAPI staining</b>		
<b>Cryptosporidium % +ve</b>	<b>100</b>	<b>%</b>
<b>Giardia % +ve</b>	<b>100</b>	<b>%</b>

**Stock specifics**

**Organism:** *Cryptosporidium parvum*  
**Strain:** Iowa  
**Source:** Bovine  
**Shed date:** 23 September 2014  
**Purification method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Organism:** *Giardia lamblia*  
**Strain:** H3  
**Source:** Gerbil  
**Shed date:** 22 October 2014  
**Purification method:** Sucrose and Percoll density gradient centrifugation

**Quality Officer:**



.....  
 Linh Le  
 25 November 2014



Accredited for compliance with ISO Guide 34.  
 Accredited Reference Material Producer  
 Accreditation Number: 14993  
 Site Number: 15276



**Water Industry Enquires:**

BTF Pty Ltd.  
 Sydney Australia  
 Tel. 61 (0) 2 8877 9150  
 Fax. 61 (0) 2 8877 9101  
[www.bffbio.com](http://www.bffbio.com)

