

**Batch Number: B 221**

**General**

**Safety:** This product is non-hazardous.  
**Storage:** 2-8°C. Do not freeze.  
**Catalogue Number:** CSCG100  
**Preparation Date:** 28 July 2016  
**Expiration Date:** 25 November 2016  
**Volume:** 1.036 ml +/- 21 ul  
**Suspension media:** Buffered saline solution.  
**Sterilisation method:** Gamma irradiation.

<b>Counts</b> <small>(Method Ref: CSCG-001)</small>		
	Mean	St.Dev.
<b>Cryptosporidium count:</b>	<b>100</b>	<b>2.1</b>
	Mean	St.Dev.
<b>Giardia count:</b>	<b>99</b>	<b>2.1</b>
<b>DAPI staining</b>		
<b>Cryptosporidium % +ve</b>	<b>100</b>	<b>%</b>
<b>Giardia % +ve</b>	<b>100</b>	<b>%</b>

**Stock specifics**

**Organism:** *Cryptosporidium parvum*  
**Strain:** Iowa  
**Source:** Bovine  
**Shed date:** 24 June 2016  
**Purification method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Organism:** *Giardia lamblia*  
**Strain:** H3  
**Source:** Gerbil  
**Shed date:** 8 July 2016  
**Purification method:** Sucrose and Percoll density gradient centrifugation

**Quality Assurance Officer:**



.....  
 Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer  
 Accreditation Number: 14993  
 Site Number: 15276



**Water Industry Enquires:**

BTF Pty Ltd.  
 Sydney Australia  
 Tel. 61 (0) 2 8877 9150  
 Fax. 61 (0) 2 8877 9101  
[www.bffbio.com](http://www.bffbio.com)

