

Batch Number: B 222

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: CSCG100
Preparation Date: 25 August 2016
Expiration Date: 23 December 2016
Volume: 1.022 ml +/- 23 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: CG-011)</small>		
	Mean	SD
Cryptosporidium count:	99	2.4
	Mean	SD
Giardia count:	100	2.0
DAPI staining		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 30 June 2016
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 10 August 2016
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Officer:



.....
 Fiona Pollock
 2 September 2016



Accredited for compliance with ISO Guide 34.
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

