

Batch Number: B 225

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: CSCG100
Preparation Date: 17 November 2016
Expiration Date: 17 March 2017
Volume: 1.047 ml +/- 13 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: CG-011)</small>		
	Mean	SD
Cryptosporidium count:	99	2.1
	Mean	SD
Giardia count:	99	1.4
DAPI staining		
Cryptosporidium % +ve	100	%
Giardia % +ve	98	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 26 October 2016
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 26 October 2016
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Officer:



.....
 Fiona Pollock
 25 November 2016



Accredited for compliance with ISO Guide 34.
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

