

Batch Number: B 252

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: CSCG100
Preparation Date: 31 May 2018
Expiration Date: 28 September 2018
Volume: 1.022 ml +/- 42 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: CG-011)</small>		
	Mean	SD
Cryptosporidium count:	99	1.3
	Mean	SD
Giardia count:	99	1.7
DAPI staining		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 11 May 2018
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 11 May 2018
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Officer



.....
 Madhuli Patil
 8 June 2018



Accredited for compliance with ISO Guide 34.
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

