

Batch Number: B 260

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: CSCG100, CS2CG100, CS4CG100, CS50CG100, W-CSCG100, W-CS4CG100
Preparation Date: 13 December 2018
Expiration Date: 12 April 2019
Volume: 1.140 ml +/- 33 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: CG-011)</small>		
	Mean	SD
Cryptosporidium count:	101	2.2
	Mean	SD
Giardia count:	101	2.4
DAPI staining		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 3 November 2018
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 23 November 2018
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Specialist



Linh Le
 2 January 2019



Accredited for compliance with ISO 17034
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:
 BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8221 9589
www.btfbio.com

