

Batch Number: B 261

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: CSCG100, CS2CG100, CS4CG100, CS50CG100, W-CSCG100, W-CS4CG100
Preparation Date: 10 January 2019
Expiration Date: 10 May 2019
Volume: 1.107 ml +/- 20 µl
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: CG-011)</small>		
	Mean	SD
Cryptosporidium count:	99	2.1
	Mean	SD
Giardia count:	100	2.5
DAPI staining		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 6 December 2018
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 14 December 2018
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Specialist



Linh Le
 21 January 2019



Accredited for compliance with ISO 17034
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquiries:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8221 9589
www.btfbio.com

