

Batch Number: B 263

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: CSCG100, CS2CG100, CS4CG100, CS50CG100, W-CSCG100, W-CS4CG100
Preparation Date: 7 March 2019
Expiration Date: 5 July 2019
Volume: 1.114 ml +/- 22 µl
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: CG-011)</small>		
	Mean	SD
Cryptosporidium count:	100	0.9
	Mean	SD
Giardia count:	100	1.9
DAPI staining		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 9 January 2019
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 19 February 2019
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Specialist



Linh Le
 15 March 2019



Accredited for compliance with ISO 17034
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquiries:
 BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8221 9589
www.bffbio.com

