

Batch Number: B 101

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: EasyPC
Preparation Date: 23 April 2014
Expiration Date: 21 August 2014
Volume: 1.996 ml +/- 17 ul
Sterilisation method: Gamma irradiation.

Counts:
(Method Ref: EPC-001)
Cryptosporidium count

Mean	St.Dev.
6406	27.9


DAPI staining
Cryptosporidium % +ve **96 %**

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Crypto Excretion Date: 11 April 2014
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

EasyPC™ is not approved for seeding water samples.

Quality Manager:



 Sue Warner



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

