

Batch Number: B 105

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: EasyPC
Preparation Date: 2 December 2014
Expiration Date: 1 April 2015
Volume: 1.989 ml +/- 13 ul
Sterilisation method: Gamma irradiation.

Counts:
(Method Ref: EPC-001)
Cryptosporidium count

Mean	St.Dev.
6362	12.3

DAPI staining
Cryptosporidium % +ve **100 %**

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Crypto Excretion Date: 14 November 2014
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

EasyPC™ is not approved for seeding water samples.

Quality Officer:



.....
 Linh Le



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

