

**Batch Number: B 111**

**General**

**Safety:** This product is non-hazardous.  
**Storage:** 2-8°C. Do not freeze.  
**Catalogue Number:** EasyPC  
**Preparation Date:** 14 July 2015  
**Expiration Date:** 11 November 2015  
**Volume:** 1.998 ml +/- 12 ul  
**Sterilisation method:** Gamma irradiation.

**Counts:**  
(Method Ref: EPC-001)  
**Cryptosporidium count**

Mean	St.Dev.
<b>6380</b>	<b>21.3</b>

**DAPI staining**  
**Cryptosporidium % +ve** **100 %**

**Stock specifics**

**Organism:** *Cryptosporidium parvum*  
**Strain:** Iowa  
**Source:** Bovine  
**Crypto Excretion Date:** 2 July 2015  
**Purification method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**EasyPC™ is not approved for seeding water samples.**

**Quality Assurance Officer:**



.....  
 Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer  
 Accreditation Number: 14993  
 Site Number: 15276



**Water Industry Enquires:**

BTF Pty Ltd.  
 Sydney Australia  
 Tel. 61 (0) 2 8877 9150  
 Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

