

Batch Number: B 113

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: EasyPC
Preparation Date: 6 October 2015
Expiration Date: 3 February 2016
Volume: 1.982 ml +/- 19 ul
Sterilisation method: Gamma irradiation.

Counts:
(Method Ref: EPC-001)
Cryptosporidium count

Mean	St.Dev.
6400	8.0

DAPI staining
Cryptosporidium % +ve **100 %**

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Crypto Excretion Date: 7 September 2015
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

EasyPC™ is not approved for seeding water samples.

Quality Assurance Officer:



Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

