

Batch Number: B 117**General**

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: EasyPC
Preparation Date: 2 May 2016
Expiration Date: 30 August 2016
Volume: 1.991 ml +/- 15 ul
Sterilisation method: Gamma irradiation.

Counts:*(Method Ref: EPC-001)***Cryptosporidium count**

Mean St.Dev.

6355 **16.4****DAPI staining****Cryptosporidium % +ve** **100 %****Stock specifics**

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Crypto Excretion Date: 13 April 2016
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

EasyPC™ is not approved for seeding water samples.**Quality Officer:**

Linh Le



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
Accreditation Number: 14993
Site Number: 15276

**Water Industry Enquires:**

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

