

Batch Number: B 120

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: EasyPC
Preparation Date: 12 September 2016
Expiration Date: 10 January 2017
Volume: 1.989 ml +/- 36 ul
Sterilisation method: Gamma irradiation.

Counts:
(Method Ref: EPC-001)
Cryptosporidium count

Mean	St.Dev.
6417	23.2

DAPI staining
Cryptosporidium % +ve **96 %**

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Crypto Excretion Date: 31 August 2016
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

EasyPC™ is not approved for seeding water samples.

Quality Assurance Officer:



.....
 Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

