

EasySeed™

- EasySeed™ is a quality control reference sample containing 100 *Cryptosporidium* and 100 *Giardia* (C&G) per tube
- Each kit contains 10 EasySeed™ tubes
- *Cryptosporidium* and *Giardia* in EasySeed™ are **inactivated** – safe to handle
- 4 month shelf life

Why use EasySeed™

- Precise number of C&G allows accurate measurement of the test recovery
- An easy way to meet **regulatory** Quality Control procedures
- **Saves labor and time** in preparing reference samples
- Inactivated organisms make it **safe** to handle in the laboratory
- **Longer shelf life** than live organisms allows flexibility in work schedules

Approvals

- United States Environmental Protection Agency
- United Kingdom Drinking Water Inspectorate
- Japan Water Works Association



EasySeed™ packaging



Water Industry Enquires:
BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

