

Batch Number: **B 635**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESC100
Preparation Date: 24 May 2018
Expiration Date: 21 September 2018
Volume: 0.598 ml +/- 4 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: ESCG-001)</small>		
	Mean	St.Dev.
Cryptosporidium count:	99	1.7
DAPI staining:		
Cryptosporidium % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 11 May 2018
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Officer



Madhuli Patil



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

