

Batch Number: **B 508**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESCG100
Preparation Date: 24 April 2014
Expiration Date: 22 August 2014
Volume: 1.245 ml +/- 12 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.


Counts		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
Cryptosporidium count:	100	1.7
	Mean	St.Dev.
Giardia count:	100	1.6
DAPI staining:		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 11 April 2014
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 09 April 2014
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Manager:



 Sue Warner



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
 www.btfbio.com

