

Batch Number: **B 536**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESCG100
Preparation Date: 23 April 2015
Expiration Date: 21 August 2015
Volume: 1.295 ml +/- 12 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts

(Method Ref: ESCG-001)

| | Mean | St.Dev. |
|-------------------------------|------------|------------|
| Cryptosporidium count: | 100 | 1.2 |

| | Mean | St.Dev. |
|-----------------------|------------|------------|
| Giardia count: | 100 | 2.0 |

DAPI staining:

Cryptosporidium % +ve **100** %

Giardia % +ve **100** %

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 22 March 2015
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 09 April 2015
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Officer:



.....
 Linh Le



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

