

Batch Number: B 578

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESCG100
Preparation Date: 13 October 2016
Expiration Date: 10 February 2017
Volume: 1.260 ml +/- 10 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
Cryptosporidium count:	99	1.5
	Mean	St.Dev.
Giardia count:	100	2.1
DAPI staining:		
Cryptosporidium % +ve	96	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 28 September 2016
Purification method: Sucrose and Percoll density gradient centrifugation

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 27 September 2016
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Officer:



Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

