

Batch Number: **B 592**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESCG100
Preparation Date: 30 March 2017
Expiration Date: 28 July 2017
Volume: 1.259 ml +/- 19 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts

(Method Ref: ESCG-001)

	Mean	St.Dev.
Cryptosporidium count:	99	1.3

	Mean	St.Dev.
Giardia count:	100	2.1

DAPI staining:

Cryptosporidium % +ve **100** %

Giardia % +ve **100** %

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 20 February 2017
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 27 February 2017
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Officer:

.....
Linh Le



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
Accreditation Number: 14993
Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

