

Batch Number: **B 595**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESCG100
Preparation Date: 27 April 2017
Expiration Date: 25 August 2017
Volume: 1.243 ml +/- 13 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
Cryptosporidium count:	99	1.5
	Mean	St.Dev.
Giardia count:	99	1.4
DAPI staining:		
Cryptosporidium % +ve	100	%
Giardia % +ve	98	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 10 March 2017
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 10 April 2017
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Officer:

Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
Accreditation Number: 14993
Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
www.btfbio.com

