

Batch Number: **B 636**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESCG100
Preparation Date: 21 June 2018
Expiration Date: 19 October 2018
Volume: 1.257 ml +/- 27 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
Cryptosporidium count:	99	1.9
	Mean	St.Dev.
Giardia count:	99	1.4
DAPI staining:		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 07 June 2018
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 05 June 2018
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Officer

Madhuli Patil



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
Accreditation Number: 14993
Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
 www.btfbio.com

