

Batch Number: B 658

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESCG100, ESCG100-5
Preparation Date: 31 January 2019
Expiration Date: 31 May 2019
Volume: 1.222 ml +/- 62 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
Cryptosporidium count:	99	2.1
	Mean	St.Dev.
Giardia count:	102	2.4
DAPI staining:		
Cryptosporidium % +ve	100	%
Giardia % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 12 January 2019
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 18 January 2019
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Manager



.....
Peter Verwoerd
 14 February 2019



Accredited for compliance with ISO 17034
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquiries:
 BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8221 9589
www.btfbio.com

