

Batch Number: **B 662**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESCG100, ESCG100-5
Preparation Date: 28 February 2019
Expiration Date: 28 June 2019
Volume: 1.287 ml +/- 8 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts

(Method Ref: ESCG-001)

	Mean	St.Dev.
Cryptosporidium count:	99	1.8

	Mean	St.Dev.
Giardia count:	99	1.6

DAPI staining:

Cryptosporidium % +ve **100 %**

Giardia % +ve **100 %**

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 09 January 2019
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Organism: *Giardia lamblia*
Strain: H3
Source: Gerbil
Shed date: 19 February 2019
Purification method: Sucrose and Percoll density gradient centrifugation

Quality Assurance Specialist



Linh Le
 15 March 2019



Accredited for compliance with ISO 17034
 Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquiries:
 BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8221 9589
www.btfbio.com

