

Batch Number: B 517

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: H-ESC-4M
Preparation Date: 12 August 2014
Expiration Date: 10 December 2014
Volume: 10.072ml +/- 36 ul
Suspension media: Buffered saline solution.
Viability: Inactivated - Gamma Irradiation

Counts

	Mean	St.Dev.
Cryptosporidium count:	3950850	98970.1

DAPI staining:

Cryptosporidium % +ve	100	%
------------------------------	------------	----------

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 04 July 2014
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Manager:



Sue Warner



Water Industry Enquires:

BTF Pty Ltd.
Sydney Australia
Tel. 61 (0) 2 8877 9150
Fax. 61 (0) 2 8877 9101
www.btfbio.com

