

Batch Number: **B 541**

**General**

**Safety:** This product is non-hazardous.  
**Storage:** 2-8°C. Do not freeze.  
**Catalogue Number:** ESC100  
**Preparation Date:** 18 June 2015  
**Expiration Date:** 16 October 2015  
**Volume:** 0.627 ml +/- 10 ul  
**Suspension media:** Buffered saline solution.  
**Sterilisation method:** Gamma irradiation.

<b>Counts</b>		
<i>(Method Ref: ESCG-001)</i>	Mean	St.Dev.
<b>Cryptosporidium count:</b>	<b>99</b>	<b>1.3</b>
<b>DAPI staining:</b>		
<b>Cryptosporidium % +ve</b>	<b>100</b>	<b>%</b>

**Stock specifics**

**Organism:** *Cryptosporidium parvum*  
**Strain:** Iowa  
**Source:** Bovine  
**Shed date:** 08 May 2015  
**Purification method:** Discontinuous sucrose and cesium chloride centrifugation gradients.

**Quality Officer:**



.....  
 Linh Le



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer  
 Accreditation Number: 14993  
 Site Number: 15276



**Water Industry Enquires:**

BTF Pty Ltd.  
 Sydney Australia  
 Tel. 61 (0) 2 8877 9150  
 Fax. 61 (0) 2 8877 9101  
[www.btfbio.com](http://www.btfbio.com)

