

Batch Number: **B 550**

General

Safety: This product is non-hazardous.
Storage: 2-8°C. Do not freeze.
Catalogue Number: ESC100
Preparation Date: 8 October 2015
Expiration Date: 5 February 2016
Volume: 0.592 ml +/- 15 ul
Suspension media: Buffered saline solution.
Sterilisation method: Gamma irradiation.

Counts <small>(Method Ref: ESCG-001)</small>		
	Mean	St.Dev.
Cryptosporidium count:	99	1.2
DAPI staining:		
Cryptosporidium % +ve	100	%

Stock specifics

Organism: *Cryptosporidium parvum*
Strain: Iowa
Source: Bovine
Shed date: 07 September 2015
Purification method: Discontinuous sucrose and cesium chloride centrifugation gradients.

Quality Assurance Officer:



 Fiona Pollock



This document is issued in accordance with NATA's accreditation requirements. Accredited for compliance with ISO Guide 34.

Accredited Reference Material Producer
 Accreditation Number: 14993
 Site Number: 15276



Water Industry Enquires:

BTF Pty Ltd.
 Sydney Australia
 Tel. 61 (0) 2 8877 9150
 Fax. 61 (0) 2 8877 9101
 www.btfbio.com

